



# Orbit Reader 20

***JAWS Commands***

**20<sup>th</sup> February, 2018  
Version 1.0.1**

# 1 JAWS Commands

## 1.1 Navigation Commands

Function	Key(s)
Move display to the left	Left
Move display to the right	Right
Move display up one line	Up
Move display down one line	Down
Pan left one braille window	Pan Left
Pan right one braille window	Pan Right
Top of active window	Select + Dots 1 2 3
Bottom of active window	Select + Dots 4 5 6
Route braille to active cursor	Dots 1 2 3 5 7 or Dots 1 2 3 5 7 Chord
Start braille auto advance mode	Select
Decrease braille auto advance speed	Select + Dot 3
Increase Braille auto advance speed	Select + Dot 6
Say previous line	Dots 2 3 or Dots 2 3 Chord
Say next line	Dots 5 6 or Dots 5 6 Chord
Say current line	Dots 1 4 or Dots 1 4 Chord
Say current word	Dots 2 5 or Dots 2 5 Chord
Say current character	Dots 3 6 or Dots 3 6 Chord
Previous document window	Dots 1 3 or Dots 1 3 Chord
Next document window	Dots 4 6 or Dots 4 6 Chord
Beginning of file	Dots 1 2 3 or Dots 1 2 3 Chord
End of file	Dots 4 5 6 or Dots 4 5 6 Chord
Page Up	Dots 1 2 3 4 5 or Dots 1 2 3 4 5 Chord
Page Down	Dots 1 2 4 5 6 or Dots 1 2 4 5 6 Chord
Home	Dots 1 2 3 4 Chord
End	Dots 1 4 5 6 Chord

Up Arrow	Dot 1 or Dot 1 Chord
Down Arrow	Dot 4 or Dot 4 Chord
Previous word	Dot 2 or Dot 2 Chord
Next word	Dot 5 or Dot 5 Chord
Previous character	Dot 3 or Dot 3 Chord
Next character	Dot 6 or Dot 6 Chord

## 1.2 Editing Commands

Function	Key(s)
Backspace (acts as shift key when typing braille characters)	Dot 7
Enter	Dot 8
Tab	Dots 4 5 or Dots 4 5 Chord
Shift + Tab	Dots 1 2 or Dots 1 2 Chord
Delete current character	Dots 1 4 5 or Dots 1 4 5 Chord
Select all	Dots 1 8 or Dots 1 8 Chord
Select from top of document	Dots 1 2 3 7 8 or Dots 1 2 3 7 8 Chord
Select to bottom of document	Dots 4 5 6 7 8
Select previous page	Dots 1 2 3 4 5 7 8
Select next page	Dots 1 2 4 5 6 7 8
Select to beginning of line	Dots 1 2 3 4 7 8
Select to end of line	Dots 1 4 5 6 7 8 or Dots 1 4 5 6 7 8 Chord
Select previous line	Dots 1 7 8
Select next line	Dots 4 7 8 or Dots 4 7 8 Chord
Select previous word	Dots 2 7 8 or Dots 2 7 8 Chord
Select next word	Dots 5 7 8 or Dots 5 7 8 Chord
Select previous character	Dots 3 7 8 or Dots 3 7 8 Chord
Select next character	Dots 6 7 8 or Dots 6 7 8 Chord
Copy to clipboard	Dots 1 4 8 or Dots 1 4 8 Chord
Cut to clipboard	Dots 1 3 4 6 8 or Dots 1 3 4 6 8 Chord

Paste clipboard	Dots 1 2 3 6 8 or Dots 1 2 3 6 8 Chord
Undo	Dots 1 3 5 6 8 or Dots 1 3 5 6 8 Chord

### 1.3 Windows Commands

Function	Key(s)
Start menu	Dots 1 3 4 7 8
Toggle menu bar	Dots 1 3 4 or Dots 1 3 4 Chord
Escape	Dots 1 3 5 6 or Dots 1 3 5 6 Chord
Alt Tab	Dots 2 3 4 5 or Dots 2 3 4 5 Chord
Minimize all applications	Dots 1 4 5 7 8

### 1.4 JAWS Commands

Function	Key(s)
Toggle keyboard help	Dots 1 4 5 6
Toggle type keys mode	Dots 1 2 3 4 5 6 Chord
JAWS cursor	Select + Dots 2 4 5
PC cursor	Select + Dots 1 2 3 4
Route JAWS cursor to PC cursor	Select + Dots 1 4
Say all to bottom	Dots 1 2 3 4 5 6 7 8
Toggle grade 2 translation	Dots 1 2 4 5 or Dots 1 2 4 5 Chord
Toggle grade 2 translation of current word	Select + Dots 2 4 5 6
Display last flash message	Select + Dots 1 2 4
JAWS find	Dots 1 2 4 Chord
JAWS find next	Dots 2 5 7 8
JAWS find previous	Dots 2 5 7 8 Chord
Toggle active cursor follows braille	Select + Dot 1
Toggle braille follows active cursor	Select + Dot 4
Cycle through braille modes	Select + Dots 1 3 4

Display 6 dot braille	Select + Dots 2 3 5
Display 8 dot braille	Select + Dots 2 3 6
Toggle cursor shape	Select + Dots 1 4 6
Restrict braille cursor	Select + Dots 1 2 3 5
Toggle characters and attributes	Select + Dots 1 6
Braille color marking	Select + Dots 3 4
Select attributes to be displayed	Select + Dots 1 3
Set 8 characters per space	Dots 2 3 6 Chord
Set unlimited characters per space	Dots 1 3 6 Chord
Toggle 8/unlimited characters per space	Select + Dots 1 3 4 6
Cycle table reading options	Select + Dots 2 3 4 5
Cycle table header options	Select + Dots 1 2 5
Graphics labeler	Select + Dots 1 2 4 5
Adjust braille options	Select + Dots 1 2
Adjust JAWS options	Select + Dots 1 2 3 6
JAWS window	Dots 2 4 5 or Dots 2 4 5 Chord
Run JAWS manager	Dots 2 3 7 8 or Dots 2 3 7 8 Chord
List of links on page	Dots 2 3 5 6 7 8
Show current time	Dots 2 3 4 5 7 8
Show script file name	Select + Dots 1 2 3 4 5
Quit JAWS	Dots 2 5 6 7 8 or Dots 2 5 6 7 8 Chord