



Orbit Reader 40

JAWS Commands

**29th December 2021
Version 1.0.1**

1 JAWS Commands

1.1 Navigation Commands

Function	Keys
Move display to the left	Left
Move display to the right	Right
Move display up one line	Up
Move display down one line	Down
Pan left one braille window	Pan Left
Pan right one braille window	Pan Right
Top of active window	Select + Dots 1 2 3
Bottom of active window	Select + Dots 4 5 6
Route braille to active cursor	Dots 1 2 3 5 7 or Dots 1 2 3 5 7 Chord
Start braille auto advance mode	Select
Decrease braille auto advance speed	Select + Dot 3
Increase Braille auto advance speed	Select + Dot 6
Say previous line	Dots 2 3 or Dots 2 3 Chord
Say next line	Dots 5 6 or Dots 5 6 Chord
Say current line	Dots 1 4 or Dots 1 4 Chord
Say current word	Dots 2 5 or Dots 2 5 Chord
Say current character	Dots 3 6 or Dots 3 6 Chord
Previous document window	Dots 1 3 or Dots 1 3 Chord
Next document window	Dots 4 6 or Dots 4 6 Chord
Beginning of file	Dots 1 2 3 or Dots 1 2 3 Chord
End of file	Dots 4 5 6 or Dots 4 5 6 Chord
Page Up	Dots 1 2 3 4 5 or Dots 1 2 3 4 5 Chord
Page Down	Dots 1 2 4 5 6 or Dots 1 2 4 5 6 Chord
Home	Dots 1 2 3 4 Chord
End	Dots 1 4 5 6 Chord
Up Arrow	Dot 1 or Dot 1 Chord
Down Arrow	Dot 4 or Dot 4 Chord
Previous word	Dot 2 or Dot 2 Chord
Next word	Dot 5 or Dot 5 Chord
Previous character	Dot 3 or Dot 3 Chord
Next character	Dot 6 or Dot 6 Chord

1.2 Editing Commands

Function	Keys
Backspace (acts as shift key when typing braille characters)	Dot 7
Enter	Dot 8
Tab	Dots 4 5 or Dots 4 5 Chord
Shift + Tab	Dots 1 2 or Dots 1 2 Chord
Delete current character	Dots 1 4 5 or Dots 1 4 5 Chord
Select all	Dots 1 8 or Dots 1 8 Chord
Select from top of document	Dots 1 2 3 7 8 or Dots 1 2 3 7 8 Chord
Select to bottom of document	Dots 4 5 6 7 8
Select previous page	Dots 1 2 3 4 5 7 8
Select next page	Dots 1 2 4 5 6 7 8
Select to beginning of line	Dots 1 2 3 4 7 8
Select to end of line	Dots 1 4 5 6 7 8 or Dots 1 4 5 6 7 8 Chord
Select previous line	Dots 1 7 8
Select next line	Dots 4 7 8
Select previous word	Dots 2 7 8
Select next word	Dots 5 7 8
Select previous character	Dots 3 7 8
Select next character	Dots 6 7 8
Copy to clipboard	Dots 1 4 8 or Dots 1 4 8 Chord
Cut to clipboard	Dots 1 3 4 6 8 or Dots 1 3 4 6 8 Chord
Paste clipboard	Dots 1 2 3 6 8 or Dots 1 2 3 6 8 Chord
Undo	Dots 1 3 5 6 8 or Dots 1 3 5 6 8 Chord

1.3 Windows Commands

Function	Keys
Start menu	Dots 1 3 4 7 8
Toggle menu bar	Dots 1 3 4 or Dots 1 3 4 Chord
Escape	Dots 1 3 5 6 or Dots 1 3 5 6 Chord
Alt Tab	Dots 2 3 4 5 or Dots 2 3 4 5 Chord
Minimize all applications	Dots 1 4 5 7 8

1.4 JAWS Commands

Function	Keys
----------	------

Toggle keyboard help	Dots 1 4 5 6
Toggle type keys mode	Dots 1 2 3 4 5 6 Chord
JAWS cursor	Select + Dots 2 4 5
PC cursor	Select + Dots 1 2 3 4
Route JAWS cursor to PC cursor	Select + Dots 1 4
Say all to bottom	Dots 1 2 3 4 5 6 7 8
Toggle grade 2 translation	Dots 1 2 4 5 or Dots 1 2 4 5 Chord
Toggle grade 2 translation of current word	Select + Dots 2 4 5 6
Display last flash message	Select + Dots 1 2 4
JAWS find	Dots 1 2 4 Chord
JAWS find next	Dots 2 5 7 8
JAWS find previous	Select + Dots 2 5 7 8
Toggle active cursor follows braille	Select + Dot 1
Toggle braille follows active cursor	Select + Dot 4
Cycle through braille modes	Select + Dots 1 3 4
Display 6 dot braille	Select + Dots 2 3 5
Display 8 dot braille	Select + Dots 2 3 6
Toggle cursor shape	Select + Dots 1 4 6
Restrict braille cursor	Select + Dots 1 2 3 5
Toggle characters and attributes	Select + Dots 1 6
Braille color marking	Select + Dots 3 4
Select attributes to be displayed	Select + Dots 1 3
Set 8 characters per space	Dots 2 3 6 Chord
Set unlimited characters per space	Dots 1 3 6 Chord
Toggle 8/unlimited characters per space	Select + Dots 1 3 4 6
Cycle table reading options	Select + Dots 2 3 4 5
Cycle table header options	Select + Dots 1 2 5
Graphics labeler	Select + Dots 1 2 4 5
Quick Settings	Select + Dots 1 2
Adjust JAWS options	Select + Dots 1 2 3 6
JAWS window	Dots 2 4 5 or Dots 2 4 5 Chord
Run JAWS manager	Dots 2 3 7 8
List of links on page	Dots 2 3 5 6 7 8
Show current time	Dots 2 3 4 5 7 8
Show script file name	Select + Dots 1 2 3 4 5
Quit JAWS	Dots 2 5 6 7 8

1.5 Cursor Routing Key Commands

Function	Keys
Move PC cursor to any cell	Any cursor router
Select block	Dot 3 with any routing key

1.6 Modifier keys and Modifier key commands

Use the modifier keys to simulate key combinations (keystrokes) that include **CTRL**, **ALT**, **WINDOWS Key**, **SHIFT**, or **the INSERT**.

To use these modifiers when typing a keystroke, do the following:

1. Hold down DOT 8 CHORD and then press the corresponding modifier keys in the keystroke. Modifier keystrokes are listed in the following table.
2. After you create the modifier portion of the keystroke, release the keys and press the remaining part of the keystroke. For example, to perform the keystroke CTRL+SHIFT+V, you press DOTS 3-7-8 CHORD, release the keys, and then press V (DOTS 1-2-3-6).

Modifier Key Name	Keystroke
INSERT	DOT 2
CTRL	DOT 3
WINDOWS Key	DOT 4
Function Keys (F1 through F12)	DOT 5
ALT	DOT 6
SHIFT	DOT 7

1.6.1 Functions Keys

To simulate the function keys (**F1** through **F12**), press **DOTS 5-8 CHORD** followed by **A** through **L** (which correspond to 1 through 12). For example, to simulate the **F6** key, press **DOTS 5-8 CHORD**, then **F (DOTS 1-2-4)**. If the function key is part of a key combination, add the appropriate modifiers while holding down **DOTS 5-8 CHORD**. For example, to simulate the keystroke **INSERT+F2**, press **DOTS 2-5-8 CHORD**, then **B (DOTS 1-2)**.

Note: Please note that Modifier keys activation shall be canceled automatically after 5 seconds if the remaining keystroke is not pressed.